

# Spill the Tea, Catechins Identified with the Hamilton PRP-C18

**Author:**

Adam L. Moore, PhD.,  
Hamilton Company



# Application

Catechins found in tea (green, white, black), are best known for their antioxidant properties. Catechins have shown to exhibit antimicrobial activity against staphylococcus, E. Coli, and H. pylori<sup>1</sup> in addition to preventing oxidative degradation. Due to these attributes, catechins have found utility in multiple arenas: commercial teas, food preservation,<sup>2</sup> free radical scavenging<sup>3</sup> and as disinfection agents.<sup>4</sup> However, some quantification protocols call for HPLC run times over 20 minutes.

In an effort to optimize the HPLC resolution after extraction of catechins from a commercial green tea bag, we have chosen a 10 mM (pH 2.5) sodium phosphate buffer for the aqueous mobile phase while utilizing acetonitrile as the organic eluent. The temperature of the separation was found to be optimized at 34°C for the best resolution and peak shape. All eight compounds

were separated in under 12 minutes. The caffeine peak, #3, is highly sensitive to temperature changes under these conditions. Increased temperatures caused increased retention, whereas cooling the column to under 20°C was observed to reduce retention to under 4 minutes while leaving the remaining peak retentions relatively unchanged. Reproducibility studies are shown of the analysis from the first injection to the 500th injection, highlighting the durability found in all Hamilton Company HPLC columns.

1. Jeon J., Kim J. H., Lee C. K. *Annals of Dermatology*. 2014, 26(5), 564.
2. He, Y. H., Shahidi, F. J. *Ag. Food Chem.* 1997, 45(11), 4262.
3. Gupta D. A., Bhaskar D. J., Gupta R. K. *Bio. Sci. Pharm. Res.* 2014, 2, 8.
4. Reysgaert, W. C. *Bio. Med. Res. Int.*, 2018 9105261.

Want more information about this and other applications?  
**Contact Hamilton today!**



# Chromatogram and Compound Results

## Column Information


Packing Material	PRP-C18 (5 µm)
Dimensions	150 x 4.6 mm
P/N	79676

## Chromatographic Conditions

Gradient	0.0–15.00 min. 10–25 %B
Temperature	34°C
Injection Volume	5 µL
Detection	UV at 214 nm
Eluent A	10 mM NaH <sub>2</sub> PO <sub>4</sub> (pH= 2.5)
Eluent B	Acetonitrile
Flow Rate	2.0 mL/min

## Compounds:

- 1: Gallocatechin
- 2: Catechin
- 3: Caffeine
- 4: Epicatechin
- 5: Epigallocatechin 3-gallate
- 6: Gallocatechin 3-gallate
- 7: Catechin 3-gallate
- 8: Epicatechin 3-gallate



[EXPLORE HAMILTON HPLC APPLICATION INDEX](#)





## About Hamilton

Hamilton Company is a global manufacturer and supplier of world-class analytical components, medical instrumentation, temperature control systems, laboratory robotics and automated liquid handling equipment. For more than 35 years, Hamilton Company has developed and manufactured pressure-stable, polymeric polystyrene-divinylbenzene (PS-DVB) HPLC columns that are used in most of the world's top chromatography labs. With a wide range of particle sizes, pore sizes, pH stability from 1 to 14, temperature resistance over 100°C, and chemistries to match most analyte types, Hamilton polymeric columns are the chromatographer's choice for challenging separations.

## Hamilton's Here to Support You



### Product Support

We are available to assist your lab with method development, optimization, troubleshooting, and product recommendations.

[CONTACT SUPPORT](#)



### Customer Service and Order Support

Get assistance with order status, lead times, and shipping information.

[CUSTOMER SERVICE](#)

# HAMILTON<sup>®</sup>

**Web:** [www.hamiltoncompany.com](http://www.hamiltoncompany.com)

**USA:** 800-648-5950

**Europe:** +40-356-635-055

#### Hamilton Americas & Pacific Rim

Hamilton Company Inc.  
4970 Energy Way  
Reno, Nevada 89502 USA  
Tel: +1-775-858-3000  
Fax: +1-775-856-7259  
[sales@hamiltoncompany.com](mailto:sales@hamiltoncompany.com)

#### Hamilton Europe, Asia & Africa

Hamilton Central Europe S.R.L.  
str. Hamilton no. 2-4  
307210 Giarmata, Romania  
Tel: +40-356-635-055  
[contact.lab.ro@hamilton-ce.com](mailto:contact.lab.ro@hamilton-ce.com)

To find a representative in your area, please visit [hamiltoncompany.com/contacts](http://hamiltoncompany.com/contacts).